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AGREED

State Security Service food products and consumer protection (conclusion of the state sanitary and epidemiological examination No. 12.2-18-2 /4369 dated 04/28/2023).

REGISTERED

Ministry of Economy of Ukraine No. 04725958/001083 dated 03/28/2023 Technical conditions of Ukraine TU U 32.5-42574330-003:2023 DCPP code 32.50.50-30.00

INSTRUCTION

for the use of an adsorbent powder of plant origin "PLANTOR" medical device "PLANTOR"

DCPP code 32.50.50-30.00, UKND code 11.160, UKT ZED code 3824 99 58 00

Product name: Adsorbent powder of plant origin "PLANTOR"

Producer: LLC "SOKOLOVSKIY-FORTUNA", Ukraine, 49100, Dnipro, Mandrykivska st. 222/21, sterile according to the recipe and technology of "SOKOLOVSKIY-FORTUNA" LLC, EDRPOU code 42574330.

Address of production facilities: Ukraine, 49100, Dnipro, Mandrykivska st., 222/21.

Release form: Powder for external use. Meets the requirements of the technical conditions of Ukraine TU U 32.5-42574330-003:2023 "Adsorbing agent in powder form of plant origin "PLANTOR". Specifications".

Pharmacological properties. Pharmacodynamics: Adsorbing, drying, antiseptic, anti-inflammatory agent for external use of plant origin. It has high hygroscopicity, adsorbs exudate, reduces vascular permeability, has a hemostatic, anti-allergic, anti-inflammatory effect, and also reduces exudation. Stimulates local immune reactions in the cell by increasing the absorbing and digesting activity of phagocytes, strengthens regenerative processes and reparative activity in the wound. A mixture of biostimulated plant extracts of juniper, bog and ginseng with essential oils of tea tree, lavender, fir, pine, thuja, sage, cedar, myrtle, eucalyptus has a pronounced antibacterial, antifungal, antimicrobial, regenerating, wound-healing effect. The product reduces the severity of local inflammation and irritation.

Pharmacokinetics: The components of the product do not penetrate through the skin and do not show a systemic effect.

Ingredients: natural polysaccharides, calcium gluconate, cationic antiseptics, plant extracts biophysically and biochemically stimulated from biologically active raw materials of plant origin with the addition of essential natural oils.

Indications for use: Stopping external bleeding of varying intensity, including damage to large venous and arterial vessels.

Recommendations for consumption: The product promotes revitalization, local coagulation of treated surfaces, stopping external bleeding of varying intensity, including damage to large venous and arterial vessels.

Contraindications: Hypersensitivity to the components of the product.

Features of application: Do not allow the product to get into the eyes.

Precautions for consumption: Individual intolerance.

Side effect: Not installed.

It is not a medicine.

Method of application and dosage.

- 1) Tear the bag and pour the powder into the wound;
- 2) On the wound with powderapply a pressure bandage;
- 3) In case of damage to large vessels of the limbs, additionally apply a tourniquet above the wound, carry out transport immobilization and evacuate the victim;
- 4) If necessary (in the case of a large wound), use several bags.

Release form: Hermetic bags of 5-20 g, it is possible, if necessary, to be equipped with a special device for introducing the agent into the wound.

Storage conditions: They are stored in transport containers in closed warehouses, at temperatures from +5 °C to +25 °C and relative air humidity (50-80%) at a distance of 1 m from heating devices.

Keep out of the reach of children.

Conditions of leave: No prescription.

Shelf Life: 60 months.

The production date and (or) expiration date, batch number are indicated on the package.